

# شركة نطاق الإنجاز للصناعات الطبية ACHIEVEMENT SCOPE COMPANY FOR MEDICAL INDUSTRIES



# ASC Medical Disinfectant Catalogue





# **ASC-SURF**

SURFACE DISINFECTANT AND CLEANER SPRAY

### أ**ي أس سي –سيرف** بخاذ مطهر ومنظف للأسطم

### **Product Description** -

ASC-SURF Spray Disinfectant is Alcohol-Free Surface Disinfectant and cleaner has been specially formulated for use in the Medical and Dental environment on non-invasive Medical devices. It is safe to use on surfaces sensitive to alcohol and contains detergent for added cleaning power and effective against wide of microorganisms.

### Product Feature -

- 100% Free from Chlorine, Phenol and Aldehyde.
- Excellent material compatibility. Materials which have been classified as safe include; Stainless Steel, Mild Steel, Aluminium, Copper, Silicon Rubber & HDPE.
- PHMB Free Modern formulation uses zero polyhexamethylene biguanide hydrochloride (PHMB).
- No-color, No-odour, Non-irritant, Non-toxic, Leaves No harmful residu & No Damaging to surface.

### **Product Composition**

Didecyl dimethylammonium chloride (DDAC),N,N-Didecyl-N-methyl-poly(oxyethyl) ammonium propionate, alkyldimethylbenzylammonium chloride, Corrosion inhibitor & Surfactant.

### Product Spectrum

- Bactericidal: Pseudomonas aeruginosa, Staphyloccocus aureus, Enterococcus hirae, Salmonella typhimurium, Legionella pneumophila, E-coli, MRSA, VRE.
- Yeasticidal: Candida albicans
- Enveloped Viruses: for HIV, Hepatitis B virus, Hepatitis C virus, Human T Cell Leukemia virus, Measles virus, Hepatitis delta virus (HDV), Influenza virus, Coronavirus (SARS, MERS), Herpesviridae, Filoviridae (Ebola), Rubella virus, Rabies virus, Flavivirus, Poxviridae.
- Tuberculocidal: Mycobacterium avium.
- Mycobactericidal: Mycobacterium terrae.

750 ml ASC-SDS-

ASC-SDS-032 -750 ml-RTU

# ASC-SURF SURFACE DISINFECTANT & CLEANER SPRAY Alcohol Free Fight teacher Alcohol Free

# **ASC-SURF**

SURFACE DISINFECTANT AND CLEANER SPRAY

### **Product Description**

ASC-SURF Spray Disinfectant is Alcohol-Free Surface Disinfectant and cleaner has been specially formulated for use in the Medical and Dental environment on non-invasive Medical devices. It is safe to use on surfaces sensitive to alcohol and contains detergent for added cleaning power and effective against wide of microorganisms.

### **Product Feature**

- 100% Free from Chlorine, Phenol and Aldehyde.
- Excellent material compatibility. Materials which have been classified as safe include; Stainless Steel, Mild Steel, Aluminium, Copper, Silicon Rubber & HDPE.
- PHMB Free Modern formulation uses zero polyhexamethylene biguanide hydrochloride (PHMB).
- No-color, No-odour, Non-irritant, Non-toxic, Leaves No harmful residu & No Damaging to surface.

### **Product Composition**

Didecyl dimethylammonium chloride (DDAC),N,N-Didecyl-N-methyl-poly(oxyethyl) ammonium propionate, alkyldimethylbenzylammonium chloride, Corrosion inhibitor & Surfactant.

### Product Spectrum

- Bactericidal: Pseudomonas aeruginosa, Staphyloccocus aureus, Enterococcus hirae, Salmonella typhimurium, Legionella pneumophila, E-coli, MRSA, VRE.
- Yeasticidal: Candida albicans
- Enveloped Viruses: for HIV, Hepatitis B virus, Hepatitis C virus, Human T Cell Leukemia virus, Measles virus, Hepatitis delta virus (HDV), Influenza virus, Coronavirus (SARS, MERS), Herpesviridae, Filoviridae (Ebola), Rubella virus, Rabies virus, Flavivirus, Poxviridae.
- Tuberculocidal: Mycobacterium avium.
- Mycobactericidal: Mycobacterium terrae.

**أي أس سـي –سـيرف** بخاخ مطمر ومنظفه ل<u>ل</u>أسطح





# **ASC-SURF**

SURFACE DISINFECTANT AND CLEANER SPRAY

أ**ي أس سي -سيرف** بخاخ مطهر ومنظف للأسطد

### **Product Description** -

ASC-SURF Spray Disinfectant is Alcohol-Free Surface Disinfectant and cleaner has been specially formulated for use in the Medical and Dental environment on non-invasive Medical devices. It is safe to use on surfaces sensitive to alcohol and contains detergent for added cleaning power and effective against wide of microorganisms.

### Product Feature -

- 100% Free from Chlorine, Phenol and Aldehyde.
- Excellent material compatibility. Materials which have been classified as safe include; Stainless Steel, Mild Steel, Aluminium, Copper, Silicon Rubber & HDPE.
- PHMB Free Modern formulation uses zero polyhexamethylene biguanide hydrochloride (PHMB).
- No-color, No-odour, Non-irritant, Non-toxic, Leaves No harmful residu & No Damaging to surface.

### **Product Composition**

Didecyl dimethylammonium chloride (DDAC),N,N-Didecyl-N-methyl-poly(oxyethyl) ammonium propionate, alkyldimethylbenzylammonium chloride, Corrosion inhibitor & Surfactant.

### Product Spectrum

- Bactericidal: Pseudomonas aeruginosa, Staphyloccocus aureus, Enterococcus hirae, Salmonella typhimurium, Legionella pneumophila, E-coli, MRSA, VRE.
- Yeasticidal: Candida albicans
- Enveloped Viruses: for HIV, Hepatitis B virus, Hepatitis C virus, Human T Cell Leukemia virus, Measles virus, Hepatitis delta virus (HDV), Influenza virus, Coronavirus (SARS, MERS), Herpesviridae, Filoviridae (Ebola), Rubella virus, Rabies virus, Flavivirus, Poxviridae.
- Tuberculocidal: Mycobacterium avium.
- Mycobactericidal: Mycobacterium terrae.

500 ml

ASC-SDS-032 -500 ml

# 

**أي أس سي –سيرف** بخاخ مطهر ومنظفه للأسطم

# **ASC-SURF**

SURFACE DISINFECTANT AND CLEANER SPRAY

### **Product Description**

ASC-SURF Spray Disinfectant is Alcohol-Free Surface Disinfectant and cleaner has been specially formulated for use in the Medical and Dental environment on non-invasive Medical devices. It is safe to use on surfaces sensitive to alcohol and contains detergent for added cleaning power and effective against wide of microorganisms.

### **Product Feature**

- 100% Free from Chlorine, Phenol and Aldehyde.
- Excellent material compatibility. Materials which have been classified as safe include; Stainless Steel, Mild Steel, Aluminium, Copper, Silicon Rubber & HDPE.
- PHMB Free Modern formulation uses zero polyhexamethylene biguanide hydrochloride (PHMB).
- No-color, No-odour, Non-irritant, Non-toxic, Leaves No harmful residu & No Damaging to surface.

### **Product Composition**

Didecyl dimethylammonium chloride (DDAC),N,N-Didecyl-N-methyl-poly(oxyethyl) ammonium propionate, alkyldimethylbenzylammonium chloride, Corrosion inhibitor & Surfactant.

### Product Spectrum

- Bactericidal: Pseudomonas aeruginosa, Staphyloccocus aureus, Enterococcus hirae, Salmonella typhimurium, Legionella pneumophila, E-coli, MRSA, VRE.
- Yeasticidal: Candida albicans
- Enveloped Viruses: for HIV, Hepatitis B virus, Hepatitis C virus, Human T Cell Leukemia virus, Measles virus, Hepatitis delta virus (HDV), Influenza virus, Coronavirus (SARS, MERS), Herpesviridae, Filoviridae (Ebola), Rubella virus, Rabies virus, Flavivirus, Poxviridae.
- Tuberculocidal: Mycobacterium avium.
- Mycobactericidal: Mycobacterium terrae.



CONCENTRATED



# Surface Disinfectant and Cleaner Solution(Concentrated)

أى أس سى-سيرڤ محلول مطهر ومنظف للأسطح (مركز)

### **Product Description**

ASC-SURF Disinfectant is Alcohol-Free Surface Disinfectant and cleaner has been specially formulated for use in the Medical and Dental environment on non-invasive Medical devices. It is safe to use on surfaces sensitive to alcohol and contains detergent for added cleaning power.

### **Product Feature**

- 100% Free from Chlorine, Phenol and Aldehyde.
- Excellent material compatibility. Materials which have been classified as safe include; Stainless Steel, Mild Steel, Aluminium, Copper, Silicon Rubber & HDPE.
- PHMB Free Modern formulation uses zero polyhexamethylene biguanide hydrochloride (PHMB).
- No-color, No-odour, Non-irritant, Non-toxic, Leaves No harmful residu & No Damaaina to surface.

### **Product Composition**

Didecyl dimethylammonium chloride (DDAC), N, N-Didecyl-N-methyl-poly(oxyethyl) ammonium propionate, alkyldimethylbenzylammonium chloride, Corrosion inhibitor &

### **Effective Spectrum & Exposure Time**

- Bactericidal 5 min (5%): Pseudomonas aeruginosa, Staphyloccocus aureus, Enterococcus hirae, Salmonella typhimurium, Legionella pneumophila, E-coli, MRSA, VRE.
- Yeasticidal 5 min (5%): Candida albicans
- Enveloped Viruses 5 min (5%): HIV, Hepatitis B virus, Hepatitis C virus and 7 min for Human T Cell Leukemia virus, Measles virus, Hepatitis delta virus (HDV), Influenza virus, Coronavirus (SARS, MERS), Herpesviridae, Filoviridae (Ebola), Rubella virus, Rabies virus, Flavivirus, Poxviridae.
- Non Enveloped Viruses 7 min (5%): Adenovirus 5, Norovirus.
- Tuberculocidal 7 min (5%): Mycobacterium avium.
- Mycobactericidal 7min (5%): Mycobacterium terrae.

5 Litter

ASC-SDS-025- CON-5L

# أى أس سى –سيرف بلس

**Surface Disinfectant Spray** بخاخ مطهر للأسطح

### **Product Description**

ASC-SURF+ Surface Disinfectant Spray has been specially formulated for use in the Medical and Dental environment. Is effective against a wide spectrum of organisms. ASC-SURF+ is a ready to use, aldehyde & phenol-free, auick disinfectant based on alcohol & it is designed for all surfaces where quick disinfection is required, e.g.:

- Patient beds, bedside tables, treatment chairs, vehicles for patient transport.
- Alcohol resistant surfaces such as bathtubs, toilet seat and transport baskets.
- Alcohol resistant non-invasive medical devices such as thermometers, sphygmomanometer and venous tourniquets.

### **Product Feature**

- 100% Free from Chlorine, Phenol and Aldehyde.
- Excellent material compatibility. Materials which have been classified as safe include; Stainless Steel, Mild Steel, Aluminium, Copper, Silicon Rubber and HDPE.
- PHMB Free Modern formulation uses zero polyhexamethylene biguanide hydrochloride (PHMB).
- No-color, No-odour, Non-irritant, Non-toxic, Leaves

### **Product Composition**

Ethanol, Propan-2-ol, Didecyl dimethylammonium chloride (DDAC), alkyldimethylbenzylammonium chloride, Corrosion inhibitor & Surfactant.

### **Effective Spectrum**

- Bactericidal: Pseudomonas aeruginosa, Staphyloccocus aureus, Enterococcus hirae, Salmonella typhimurium, Legionella pneumophila, E-coli, MRSA, VRE.
- Yeasticidal: Candida albicans
- Enveloped Viruses: for HIV, Hepatitis B virus, Hepatitis C virus, Human T Cell Leukemia virus, Measles virus, Hepatitis delta virus (HDV), Influenza virus, Coronavirus (SARS, MERS), Herpesviridae, Filoviridae (Ebola), Rubella virus, Rabies virus, Flavivirus, Poxviridae.
- Tuberculocidal: Mycobacterium avium.
- Mycobactericidal: Mycobacterium terrae







# اي اس سي -سيرف بلس (+) **ASC-SURF**

Surface Disinfectant Spray ىخاخ مطهر للأسطح

### Product Description

ASC-SURF+ Surface Disinfectant Spray has been specially formulated for use in the Medical and Dental environment. Is effective against a wide spectrum of organisms. ASC-SURF+ is a ready to use, aldehyde & phenol -free, quick disinfectant based on alcohol & it is designed for all surfaces where quick disinfection is required, e.g.:

- Patient beds, bedside tables, treatment chairs, vehicles for patient transport.
- Alcohol resistant surfaces such as bathtubs, toilet seat and transport baskets.
- Alcohol resistant non-invasive medical devices such as thermometers, sphygmomanometer and venous tourniquets.

### Product Feature

- 100% Free from Chlorine, Phenol and Aldehyde.
- Excellent material compatibility. Materials which have been classified as safe include; Stainless Steel, Mild Steel, Aluminium, Copper, Silicon Rubber and HDPE.
- PHMB Free Modern formulation uses zero polyhexamethylene biguanide hydrochloride (PHMB).
- No-color, No-odour, Non-irritant, Non-toxic, Leaves

### **Product Composition**

Ethanol, Propan-2-ol, Didecyl dimethylammonium chloride (DDAC), alkyldimethylbenzylammonium chloride, Corrosion inhibitor & Surfactant.

### **Effective Spectrum**

- Bactericidal: Pseudomonas aeruginosa, Staphyloccocus aureus, Enterococcus hirae, Salmonella typhimurium, Legionella pneumophila, E-coli, MRSA, VRE.
- Yeasticidal: Candida albicans
- Enveloped Viruses: for HIV, Hepatitis B virus, Hepatitis C virus, Human T Cell Leukemia virus, Measles virus, Hepatitis delta virus (HDV), Influenza virus, Coronavirus (SARS, MERS), Herpesviridae, Filoviridae (Ebola), Rubella virus, Rabies virus, Flavivirus, Poxviridae.
- Tuberculocidal: Mycobacterium avium.
- Mycobactericidal: Mycobacterium terrae.

500 ml

ASC-SDS-PLUS-500ml-RTU



Fast acting to eliminates viruses, germs, bacteria and fungi

،سريع المفعول للقضاء على الفيروسات والجراثيم البكتيريا والفطريات

### **ASC-CHLOR**

**CHLORINE DISINFECTANT TABLETS** 

### **Product Description**

Chlorine tablets from asc are considered an ideal solution for sterilization and it also dissolves quickly. chlorine tablets are used to disinfect blood spills and general environmental disinfection. chlorine tablets are stable, compact and easy to store for longer periods until required for use, without the spill, contamination and storage risks associated with liquid chlorine solutions.

### **Ingredients**

TRICHLORO-S-TRIAZINETRIONEA: .....99.05% OTHER INGREDINTS:..... 0.95%

Dilusion —				
REQUIRED CONCENTRATION OF AVAILABLE CHLORINE (PPM)	Ditusion			
	NO. OF TABLETS	0.5G TABLET/ WATER	1.7G TABLET/ WATER	2.5G TABLET/ WATER
10000	4	120ml	300ml	600ml
5000	2	120ml	300ml	600ml
2500	1	120ml	300ml	600ml
1000	1	300ml	750ml	1.5 L
500	1	600ml	1.5 L	3 L
50	1	6 L	15 L	30 L

# BLOOD& BODY FLUID SPILS 10000 PPM LABORATORY DISCARD JARS 2500 PPM STAINLESS STEEL INSTRUMENTS 600 PPM TOILETS, DRAINS AND SINKS 400 PPM FOOD PREPARATION, SURFACES AND FLOORS 200 PPM DISHCLOTHS, MOPS, ETC. 50 PPM





### **CHLOROSEPT 4%**

### Hand Disinfectant

- Chlorhexidine gluconate 4 % W/V.
   Isopropyl Alcohol 4 %. W/V

### **Product Description**

ASC-CHLOROSEPT, is hand disinfectant has been specially formulated for use in the Medical and Dental envirinments. It is a hand Disinfectant that provides 100% better protection. A gentle gel to foam formula deep cleans & disinfects without leaving skin

أى أس سى كلوروسيبت صمم خصيصا للإستخدام في البيئة الطبية وطب الأسنان. إنه مطهر لليدين يوفر حماية .أفضل بنسبة ١٠٠ ٪. تركيبة من الجل اللطيف إلى الرغوى تطهر وتنظف بعمق دون ترك البشرة جافة أو مشدودة

### Active Ingredients

- Chlorhexidine aluconate 4 % W/V.
- Isopropyl Alcohol 4 %. W/V

1000 ml

ASC-CHL-4-1000 ml

500 ml

ASC-CHL-4-500 ml



### **CHLOROSEPT 2%**

كلوروسيبت %2

### Hand Disinfectant

- Chlorhexidine gluconate 2 % W/V.
  Isopropyl Alcohol 4 %. W/V

### **Product Description**

ASC-CHLOROSEPT, is hand disinfectant has been specially formulated for use in the Medical and Dental envirinments. It is a hand Disinfectant that provides 100% better protection. A gentle gel to foam formula deep cleans & disinfects without leaving skin

أى أس سي كلوروسيبت صمم خصيصا للإستخدام في البيئة الطبية وطب الأسنان. إنه مطهر لليدين يوفر حماية .أفضل بنسبة ١٠٠ ٪. تركيبة من الجل اللطيف إلى الرغوي تطهر وتنظف بعمق دون ترك البشرة جافة أو مشدودة

### - Active Ingredients

- Chlorhexidine gluconate 2 % W/V.
- Isopropyl Alcohol 4 %. W/V

1000 ml

ASC-CHL-2-1000 ml

500 ml

ASC-CHL-2-500 ml



### **CHLOROSEPT 1.5%**

كلوروسيست 1.5%

### - Product Description

ASC-CHLOROSEPT, is hand disinfectant has been specially formulated for use in the Medical and Dental envirinments. It is a hand Disinfectant that provides 100% better protection. A gentle gel to foam formula deep cleans & disinfects without leaving skin

أي أس سي كلوروسيبت صممر خصيصا للإستخدام في البيئة الطبية وطب الأسنان. إنه مطهر لليدين يوفر حماية .أفضل بنسبة ١٠٠ ٪. تركيبة من الجل اللطيف إلى الرغوي تطهر وتنظف بعمق دون ترك البشرة جافة أو مشدودة

### **Active Ingredients**

- Chlorhexidine aluconate 1.5% W/V.
- Cetrimide 15 %. W/V

1000 ml

ASC-CHL-1.5-1000 ml

### Hand Disinfectant

- Chlorhexidine gluconate 1.5 % W/V.
  Cetrimide 15 %. W/V





### **CHLOROSEPT 4%**

کلوروسیبت %4

### Hand Disinfectant

- Chlorhexidine gluconate 4 % W/V.
  Isopropyl Alcohol 4 %. W/V

### · Product Description

ASC-CHLOROSEPT, is hand disinfectant has been specially formulated for use in the Medical and Dental envirinments. It is a hand Disinfectant that provides 100% better protection. A gentle gel to foam formula deep cleans & disinfects without leaving skin dry or tight.

أى أس سي كلوروسيبت صمم خصيصا للإستخدام في البيئة الطبية وطب الأسنان. إنه مطهر لليدين يوفر حماية .أفضل بنسبة ١٠٠ ٪. تركيبة من الجل اللطيف إلى الرغوي تطهر وتنظف بعمق دون ترك البشرة جافة أو مشدودة

### - Active Ingredients -

- Chlorhexidine aluconate 4 % W/V.
- Isopropyl Alcohol 4 %. W/V

1000 ml

ASC-CHL-4-B - 1000 ml



### **CHLOROSEPT 2%** كلوروسيبت 12

### Hand Disinfectant

- Chlorhexidine gluconate 2 % W/V.
  Isopropyl Alcohol 4 %. W/V

### **Product Description**

ASC-CHLOROSEPT, is hand disinfectant has been specially formulated for use in the Medical and Dental envirinments. It is a hand Disinfectant that provides 100% better protection. A gentle gel to foam formula deep cleans & disinfects without leaving skin dry or tight.

أى أس سي كلوروسيبت صمم خصيصا للإستخدام في البيئة الطبية وطب الأسنان. إنه مطهر لليدين يوفر حماية .أفضل بنسبة ١٠٠ ٪. تركيبة من الجل اللطيف إلى الرغوى تطهر وتنظف بعمق دون ترك البشرة جافة أو مشدودة

### Active Ingredients

- Chlorhexidine gluconate 2 % W/V.
- Isopropyl Alcohol 4 %. W/V

1000 ml

ASC-CHL-2-B - 1000 ml





### CHLOROSEPT کلوروســــيبت

**Hond** sanitizer gel

### **Product Description** -

Chlorosept HAND SANITIZER Gel kills 99.9% of common germs within 30 seconds. just rub into hands and let dry, enhanced with moisturizers specially glycerine & vitamin E to help moisturize the skin.

:وصف المنتج

كلوروسيبت جل مطهر اليدين يقتل ٪ ٩٩،٩ من الجراثيم في غضون ٣٠ ثانية. ضعه على اليدين واتركه حتي .يجف ، معزز بالمرطبات خاصة الجلسرين و فيتامين ه للمساعدة في ترطيب البشرة

### - Active Ingredients

Water(aqua), Ethyl alcohol, Vitamin E, Glycerine and Carbomer.

KILLS 99.9% OF BACTERIA &VIRUSES

600 ml

ASC-CHL-HSG-F- 600ml



### CHLOROSEPT کلوروســــيبت

Hand sanitizer gel

### **Product Description**

ASC-CHLOROSEPT, is hand disinfectant has been specially formulated for use in the Medical and Dental envirinments. It is a hand Disinfectant that provides 100% better protection. A gentle gel to foam formula deep cleans & disinfects without leaving skin

:وصف المنتج

أي أس سي كلوروسيبت صمم خصيصا للإستخدام في البيئة الطبية وطب الأسنان. إنه مطهر لليدين يوفر حماية .أفضل بنسبة ١٠٠ ٪. تركيبة من الجل اللطيف إلى الرغوي تطهر وتنظف بعمق دون ترك البشرة جافة أو مشدودة

### - Active Ingredients -

- Chlorhexidine gluconate 2 % W/V.
- Isopropyl Alcohol 4 %. W/V



1000 ml

ASC-CHL-HSG- 1000ml





# CHLOROSEPT کلوروسیبت

Hand sanitizer gel

### **Product Description** -

Chlorosept HAND SANITIZER Gel kills 99.9% of common germs within 30 seconds, just rub into hands and let dry, enhanced with moisturizers specially glycerine & vitamin E to help moisturize the skin.

:وصف المنتج

كلوروسيبت جل مطهر اليدين يقتل ٪ ٩٩،٩ من الجراثيم في غضون ٣٠ ثانية. ضعه على اليدين واتركه حتي .يجف ، معزز بالمرطبات خاصة الجلس ين و فيتامين ه للمساعدة في ترطيب البشرة

### - Active Ingredients

Water(aqua), Ethyl alcohol, Vitamin E, Glycerine and Carbomer.



KILLS 99.9% OF BACTERIA &VIRUSES

1000 ml

ASC-CHL-HSG-B- 1000ml

# CHLOROSEPT کلوروسیبت

Hand sanitizer gel

### **Product Description**

ASC-CHLOROSEPT, is hand disinfectant has been specially formulated for use in the Medical and Dental envirinments. It is a hand Disinfectant that provides 100% better protection. A gentle gel to foam formula deep cleans & disinfects without leaving skin

:وصف المنتج

أي أس سي كلوروسيبت صمم خصيصا للإستخدام في البيئة الطبية وطب الأسنان. إنه مطهر لليدين يوفر حماية .أفضل بنسبة ۱۰۰ ٪. تركيبة من الجل اللطيف إلى الرغوي تطهر وتنظف بعمق دون ترك البشرة جافة أو مشدودة

### - Active Ingredients -

- Chlorhexidine gluconate 2 % W/V.
- Isopropyl Alcohol 4 %. W/V



60 ml

ASC-CHL-HS- 60ml





Isopropyl Alcohol 70 %

### أي أس سي-أيزو کحول أيزوبروبيل ٧٠٪

### **Product Description**

ASC-ISO Isopropyl Alcohol 70% is Rubbing Alcohol and it is the perfect solution for cleaning and sanitizing skin and surfaces. The 70% isopropyl alcohol solution is a powerful agent for cleaning cuts, burns, scrapes, and scratches while also promoting a speedy recovery. Disinfectant sanitizer alcohol for every day, all-purpose cleaning. Powerful Surface Disinfectant, effective against wide range of Bacteria & Virus. Effective for Cleaning Glass & Mirrors, Evaporates Quickly.

### وصف المنتج

أي أس سي-أيزو هو الحل الأمثل لتنظيف وتعقيم البرة والأسطح يعتر محلول كحول الأيزوبروبيل %70 عاماً قويًا لتنظيف الأغراض خال اليوم tالجروح والحروق والخدوش مع تعزيز الشفاء العاجل أيضًا. محلول مطهر ومنظف لكل مطهر قوي لأسطح ، فعال ضم مجموعة واسعة من البكتيرياوالفيروسات. فعال في تنظيف النظم الليامية في سيامة الزجاج والمرايا ، يتبخر بسرعة

### - Indications

- For antiseptic purposes, disinfection of hands and surfaces.
- For the treatment of minor cuts & Abrasions.
- For Rubbing & massaging.

### - Product Composition -

- ACTIVE INGREDIENTS: ISOPROPYL A L C O H O L .
- INACTIVE INGREDIENTS: Purified water.

250 ml

ASC-ISO-029 - 250 ml



Isopropyl Alcohol 99.9 %

Fast acting to eliminates viruses, germs, bacteria and fungi

،سريع المفعول للقضاء على الفيروسات والجراثيم البكتيريا والفطريات

### Product Description

ASC-ISO Isopropyl Alcohol 99.9 % is Rubbing Alcohol and it is the perfect solution for cleaning and sanitizing skin and surfaces. The 99.9% isopropyl alcohol solution is a powerful agent for cleaning cuts, burns, scrapes, and scratches while also promoting a speedy recovery. Disinfectant sanitizer alcohol for every day, all-purpose cleaning. Powerful Surface Disinfectant, effective against wide range of Bacteria & Virus. Effective for Cleaning Glass & Mirrors, Evaporates Quickly.

### وصف المنتج

أي أس سي-أيزو هو الحل الأمثل لتنظيف وتعقيم البشرة والأسطح. يعتبر محلول كحول الأيزوبروبيل 99.9% عام الاً قو يا لتنظيف الجروح والحروق والخدوش مع تعزيز الشفاء العاجل أيضًا. محلول مظهر ومنظف .لكل الأغراض خلال اليوم . مظهر قوي للأسطح ، فعال ضد مجموعة واسعة من البكتيريا والفيروسات فعال في تنظيف الزجاج والمرايا ، يتبخر بسرعة

### - Indications

- For antiseptic purposes, disinfection of hands and surfaces.
- For the treatment of minor cuts & Abrasions.
- · For Rubbing & massaging.

### Product Composition

• ISOPROPYL ALCOHOL 99.9 %.



**5 Litter** 

ASC-ISO-030 - 5L



# ASC-HYDRO Hydrogen Peroxide 3 %

### **Product Description**

ASC-HYDRO2 Disinfectant is Alcohol-Free Disinfectant has been specially formulated for use in the Medical and Dental environment. highly effective Disinfectant act as Bactericidal, Fungicidal, Viricidal and Sporicidal. It is safe to use on surfaces sensitive to alcohol and contains detergent for added cleaning power.

### وصف المنتج

أى أس سى – هيدرو 2 مطهر صمم خ يطًّا ص للاستخدام في البيئة الطبية وطب الأسنان مادة فعالة للغاية تستخدم كمضاد للبكتريا والفطريات والفيروسات والجراثيم. المنتج خالى من الكحول و فعال ضد مجموعة وا عة ن الكائ ات الح ة الدقية

### - Indications

Hydrogen peroxide as active ingredient, Colloidal silver, purified water.

### **Product Application** -

- Disinfectant
- Sanitizer
- Bleaching Agent
- Surface Cleaner

250 ml

ASC-HPO6 - 007 - 250 ml



# ASC-HYDRO Hydrogen Peroxide 3 % STABILIZED

### **Product Description**

ASC-ISO Isopropyl Alcohol 99.9 % is Rubbing Alcohol and it is the perfect solution for cleaning and sanitizing skin and surfaces. The 99.9% isopropyl alcohol solution is a powerful agent for cleaning cuts, burns, scrapes, and scratches while also promoting a speedy recovery. Disinfectant sanitizer alcohol for every day, all-purpose cleaning. Powerful Surface Disinfectant, effective against wide range of Bacteria & Virus. Effective for Cleaning Glass & Mirrors, Evaporates Quickly.

### وصف المنتج

أي أس سي-أيزو هو الحل الأمثل لتنظيف وتعقيم البشرة والأسطح. يعتبر مجلول كحول الأيزوبروبيل 99.9% عام ا لا قو يا لتنظيف الجروح والحروق والخدوش مع تعزيز الشفاء العاجل أيضًا. محلول مطهر ومنظف .لكل الأغراض خلال اليوم . مطهر قوي للأسطح ، فعال ضد مجموعة واسعة من البكتيريا والفيروسات فعال في تنظيف الزجاج والمرايا ، يتبخر بسرعة

### - Indications

- For antiseptic purposes, disinfection of hands and surfaces.
- For the treatment of minor cuts & Abrasions.
- For Rubbing & massaging.

### **Product Composition**

• ISOPROPYL ALCOHOL 99.9 %.



300 ml

ASC-HPO6 - SP - 300 ml



# ASC-HYDRO Hydrogen Peroxide 5 %

### **Product Description**

ASC-HYDRO2 Disinfectant is Alcohol-Free Disinfectant has been specially formulated for use in the Medical and Dental environment, highly effective Disinfectant act as Bactericidal, Fungicidal, Viricidal and Sporicidal. It is safe to use on surfaces sensitive to alcohol and contains detergent for added cleaning power.

### وصف المنتج

أى أس سى – هيدرو 2 مطهر صمم خ يصًّا ص للاستخدام في البيئة الطبية وطب الأسنان مادة فعالة للغاية تستخدم كمضاد للبكتريا والفطريات والفيروسات والجراثيم. المنتج خالى من الكحول و فعال ضد مجموعة وا عة ن الكائ ات الح ة الدقية

### · Indications

Hydrogen peroxide as active ingredient, Colloidal silver, purified water.

### **Product Application** -

- Disinfectant
- Sanitizer
- Bleaching Agent
- Surface Cleaner

750 ml ASC-HPO6 - 007 - 1000 ml



# ASC-HYDRO Hydrogen Peroxide 5 % STABILIZED 2

### Product Description

ASC-ISO Isopropyl Alcohol 99.9 % is Rubbing Alcohol and it is the perfect solution for cleaning and sanitizing skin and surfaces. The 99.9% isopropyl alcohol solution is a powerful agent for cleaning cuts, burns, scrapes, and scratches while also promoting a speedy recovery. Disinfectant sanitizer alcohol for every day, all-purpose cleaning. Powerful Surface Disinfectant, effective against wide range of Bacteria & Virus. Effective for Cleaning Glass & Mirrors, Evaporates Quickly.

### وصف المنتج


أي أس سي-أيزو هو الحل الأمثل لتنظيف وتعقيم البشرة والأسطح. يعتبر مجلول كحول الأيزوبروبيل 99.9% عام ا لا قو يا لتنظيف الجروح والحروق والخدوش مع تعزيز الشفاء العاجل أيضًا. محلول مطهر ومنظف .لكل الأغراض خلال اليوم . مطهر قوي للأسطح ، فعال ضد مجموعة واسعة من البكتيريا والفيروسات فعال في تنظيف الزجاج والمرايا ، يتبخر بسرعة

### - Indications

- For antiseptic purposes, disinfection of hands and surfaces.
- For the treatment of minor cuts & Abrasions.
- For Rubbing & massaging.

### **Product Composition**

• ISOPROPYL ALCOHOL 99.9 %.



1000 ml

ASC-HPO6 - SP - 1000 ml



# ASC-HYDRO Hydrogen Peroxide 5 %

### **Product Description**

ASC-HYDRO2 Disinfectant is Alcohol-Free Disinfectant has been specially formulated for use in the Medical and Dental environment. highly effective Disinfectant act as Bactericidal, Fungicidal, Viricidal and Sporicidal. It is safe to use on surfaces sensitive to alcohol and contains detergent for added cleaning power.

### وصف المنتج

أى أس سى – هيدرو 2 مطهر صمم خ يضًّا ص للاستخدام في البيئة الطبية وطب الأسنان مادة فعالة للغاية تستخدم كمضاد للبكتريا والفطريات والفيروسات والجراثيم. المنتج خالى من الكحول و فعال ضد مجموعة وا عة ن الكائ ات الح ة الدقية

### - Indications

Hydrogen peroxide as active ingredient, Colloidal silver, purified water.

### Product Application -

- Disinfectant
- Sanitizer
- Bleaching Agent
- Surface Cleaner

5 Litter

ASC-HPO5-007 - 5 L



# ASC-HYDRO Hydrogen Peroxide 1.41 %

### - Product Description

ASC-HYDRO2 Disinfectant is Alcohol-Free Disinfectant has been specially formulated for use in the Medical and Dental environment. highly effective Disinfectant act as Bactericidal, Fungicidal, Viricidal and Sporicidal. It is safe to use on surfaces sensitive to alcohol and contains detergent for added cleaning power.

### وصف المنتج

أى أس سى – هيدرو 2 مطهر صمم خ يصًّا ص للاستخدام في البيئة الطبية وطب الأسنان مادة فعالة للغاية تستخدم كمضاد للبكتريا والفطريات والفيروسات والجراثيم. المنتج خالى من الكحول و فعال ضد مجموعة وا عة ن الكائ ات الح ة الدقية

### · Indications

Hydrogen peroxide as active ingredient, Colloidal silver, purified water.

### **Product Application**

- Disinfectant
- Sanitizer
- Bleaching Agent
- Surface Cleaner

1000 ml ASC-HPO- B - 1000 ml



### للصناعات الطبية شركة نطاق الإنجاز

Disinfectant solution for Dental and Medical Instruments (Concentrated)

أى أس سى-انسترو محلول مطهر للأدوات الطبية وأدوات طب الأسنان (مركز)

### **Product Description**

ASC-INSTRU is a brand-new Instrument Disinfectant and Cleaner concentrate solution ideally suited for the manual disinfection of instruments, burs and endodontic instruments and items prior to sterilization.

ً انسرو محلول مطهر ومنظف عالي التركيز صُمم خصيصاً لأدوات الطبية يعمل بشكل مثالي لتطهر الأدوات قبل عملية التعقيم.

### - Product Feature

- 100% Free from Chlorine, Phenol, and Aldehyde
- Excellent material compatibility. Materials which have been classified as safe include; Stainless Steel, Mild Steel, Aluminium, Copper, Silicon Rubber and HDPE.
- PHMB Free Modern formulation uses zero polyhexamethylene biguanide hydrochloride
- Non-toxic, Leaves No harmful residu & No Damaging to surface.

### **Product Composition**

Didecyl dimethylammonium chloride (DDAC), N,N-Didecyl-N-methylpoly(oxyethyl) ammonium propionate, alkyldimethylbenzylammonium chloride, Corrosion inhibitor, Tetrasodium ethylene diamine tetraacetate, Surfactant, Perfume and Dye.

### Effective Spectrum & Exposure Time

- Bactericidal 1 min (5%): Pseudomonas aeruginosa, Staphyloccocus aureus, Enterococcus hirae, Salmonella typhimurium, Legitonella pneumophila, E-coli, MRSA, VRE
- Yeasticidal 1 min (5%): Candida albicans Enveloped Viruses 5 min (5%): HIV. Hepatitis B virus. Hepatitis C virus, Human T Cell Leukemia virus, Measles virus, Hepatitis delta virus (HDV), Influenza virus, Coronavirus (SARS, MERS), Herpesviridae, Filoviridae (Ebola), Rubella virus, Rabies virus, Flavivirus, Poxviridae
- Non Enveloped Viruses 5 min (5%): Adenovirus 5, Norovirus
- Tuberculocidal 5 min (5%): Mycobacterium avium
- Mycobactericidal 5 min (5%): Mycobacterium terrae



5 Litter ASC-INSD-024- 5L 1000 ml ASC-INSD-024- 1L

Disinfectant & Cleaner for Surface & Floor

أى أس سى-فلور مطهر ومنظف للأسطح والأرضيات

### **Product Description**

ASC-FLOR floor disinfectant works as a cleaner, deodorizer and viricide in one easy step. It Fights dangerous bacteria and viruses for a cleaner and healthier workspace. Phoshate-free formula won,t damage floors. Easy-to-use formula so you can clean floors without all the fuss. Suitable for floors, walls, toilet, drains, sinks, rubbish bins, countertop and kitchen benches. Suitable for use on all water-washable surfaces such as ceramic, stainless steel & plastic.

وصف المنتج طوة واحدة سهلة. يحارب البكت ريريا والف ريروسات للحصول على مساحة عمل

سور من سبب عرب المراحية المركزين وكوركست وكوركست من المركزية المر

### **Product Feature**

- 100% Free from Chlorine, Phenol and Aldehyde.
- Excellent material compatibility. Materials which have been classified as safe include; Stainless Steel, Mild Steel, Aluminium, Copper, Silicon Rubber and HDPE.
- PHMB Free Modern formulation uses zero polyhexamethylene biguanide hydrochloride (PHMB).

### Ingredients

Didecyl dimethylammonium chloride (DDAC), N,NDidecyl-N-methyl-poly(oxyethyl) ammoniumpropionate, alkyldimethylbenzylammonium chloride, Corrosion inhibitor & Surfactant.

### Effective Spectrum -

Bactericide, fungicide, mycobactericidal and viricide.

### **Economical**

The concentrated disinfectant has an economical dilution ratio of only 8% or 10%.





### **ASC-HYPO**

5.25 Sodium Hypochlorite صودیوم هیبوکلوریت ۲٫۵۵٪

### **Product Description**

ASC-HYPO 5.25 % is a broad-spectrum disinfectant that effectively against viruses, bacteria, fungi, and mycobacterium. It is the solution for the irrigation, debridement, and disinfection of root canals during endodontic treatment. It prevents from the teeth discoloration that may occur after filling of the root canal which was not rinsed.

### وصف المنتج

أي أس سي هيبو ٥.٢0 % مطهر واسع المدى فعال ضد الفيروسات والبكتريا و الفطريات. يستخدم ك محلول للتروية و ازالة التسبخات و تطهير قنوات العصب اثناء علاج الجذور . ويمنع الاسنان من التصبغات وتغيير لونها .الذى عادة ما يحدث بعد حشو العصب اذا لم يتم استخدامه

### **Product Composition**

Contains Sodium Hypochlorite as the main active ingredient.

**5 Litter** ASC-SH3-011 - 5L

1000 ml ASC-SH5 - 013 - 1000 ml



### ASC- Rust & Lime Scale Remover

### مزيل للصدأ

Design to ease the process of rust & lime scale Removal from the surface of instruments

### **Product Description**

ASC- RUST REMOVAL is a lightly foaming acidic cleaner that will prolong the working life and improve the appearance of adected instruments designed for the cleaning, maintenance and removal of rust and corrosion from stainless steel surgical instrument.

### وصف المنتج

منتج إزالة الصدأ من )أي أس سي( منظف حمي برغوة منخفضة يطيل من عمر الأدوات الجراحية. تم تصميمه لتنظيف وصيانة و إزالة الصدأ والتآكل من الأدوات الجراحية المصنوعة من الفولاذ.

### - Product Feature

- Dissolve lime scale easily.
- Remove rust stains.
- Quick acting.
- Remove tough inorganic stains.
- Safe to use on most surfaces.hylene biguanide hydrochloride (PHMB).

### Ingredients -

Phosphoric Acid, Anti-corrosion & Auxiliary material.

ASC-Rust & Line Scale
Remover
DESCRIPTION DURFNE PROCESS OF
SUBSTANCES CONSTRUCTION
SUBSTANCES OF RESOURCES

SUBSTANCES O

5 L

ASC - RR - 070 - 5L

### للصناعات الطبية شركة نطاق الإنجاز

### ASC-JET-7 Weekly Dental Aspirator Disinfectant and Cleaner (Concentrated)

أي أس سي - چيت- 7 مطهر ومنظف مركز لأدوات طب الأسنان وأجهزة الشفط - أسبوعي

### **Product Description**

ASC-JET-7 Weekly Dental Aspirator Used for the daily cleaning and disinfecting of aspirator systems, can be used on wet line, dry line suction systems. Highly effective cleaning and disinfecting product. Low-foaming and non-abrasive aldehyde, chlorine and phenol free solution. Easy to use and measure. Effective against a wide spectrum of organisms. Alkaline based solution for removal of Weekly deposits.

### **Product Feature**

- 100% Free from Chlorine, Phenol, Aldehyde and Non-Foaming Excellent material compatibility. Materials which have been classified as safe include; Stainless Steel, Mild Steel, Aluminium, Copper, Silicon Rubber and HDPE.
- PHMB Free Modern formulation uses zero polyhexamethylene biguanide hydrochloride

### **Product Composition**

Didecyl dimethylammonium chloride (DDAC), alkyldimethylbenzylammonium chloride, Anti-limescale, Tetrasodium ethylene diamine tetraacetate, Nonionic surfactant, Antifoam, Perfume and Dve.

### Effective Spectrum & Exposure Time -

- Bactericidal 15 min (5%): Pseudomonas aeruginosa, Staphyloccocus aureus, Enterococcus hirae, Salmonella typhimurium, Legionella pneumophila, E-coli, MRSA, VRE.
- Yeasticidal 15 min (5%): Candida albicans
- Enveloped Viruses (5%): 5 min for HIV, Hepatitis B virus, Hepatitis C virus and 30 min for Human T Cell Leukemia virus, Measles virus, Hepatitis delta virus (HDV), Influenza virus, Coronavirus (SARS, MERS), Herpesviridae, Filoviridae (Ebola), Rubella virus, Rabies virus, Flavivirus, Poxviridae.

ASC-JET-7 Weekly Dental Aspirator

5 Litter

ASC-JET-7-026-D-5L

## ASC-JET Daily Dental Aspirator Disinfectant and Cleaner (Concentrated)

أى أس سى - چيت مطهر ومنظف مركز لأدوات طب الأسنان وأجهزة الشفط - يومى

### Product Description

ASC-JET Daily Dental Aspirator Used for the daily cleaning and disinfecting of aspirator systems, can be used on wet line, dry line suction systems. Highly effective cleaning and disinfecting product. Low-foaming and non-abrasive aldehyde, chlorine and phenol free solution. Easy to use and measure. Effective against a wide spectrum of organisms. Alkaline based solution for removal of daily deposits.

### **Product Feature**

- 100% Free from Chlorine, Phenol, Aldehyde and Non-Foaming Excellent material compatibility. Materials which have been classified as safe include; Stainless Steel, Mild Steel, Aluminium, Copper, Silicon Rubber and HDPE.
- PHMB Free Modern formulation uses zero polyhexamethylene biguanide hydrochloride (PHMB).

### **Product Composition**

Didecyl dimethylammonium chloride (DDAC), alkyldimethylbenzylamtmonium chloride, Anti-limescale, Tetrasodium ethylene diamine tetraacetate, Nonionic surfactant, Antifoam, Perfume and Dye.

### Effective Spectrum & Exposure Time

- Bactericidal 15 min (5%): Pseudomonas aeruginosa, Staphyloccocus aureus, Enterococcus hirae, Salmonella typhimurium, Legionella pneumophila, E-coli, MRSA, VRE.
- Yeasticidal 15 min (5%): Candida albicans.
- Enveloped Viruses (5%): 5 min for HIV, Hepatitis B virus, Hepatitis C virus and 30 min for Human T Cell Leukemia virus, Measles virus, Hepatitis delta virus (HDV), Influenza virus, Coronavirus (SARS, MERS), Herpesviridae, Filoviridae (Ebola), Rubella virus, Rabies virus, Flavivirus, Poxviridae.
- Non Enveloped Viruses 30 min (5%): Adenovirus 5, Norovirus.
- Tuberculocidal 60 min (5%): Mycobacterium avium.
- Mycobactericidal 60 min (5%): Mycobacterium terrae.





# FLÓRÓSEPT

PERFUMED - HAND SOAP -

### **Product Description**

Liquid Florosept-handwash is a ready to use, detergent-based hand wash that provides superior cleaning and anti-bacterial effects. Florosept-handwash is strong and very effectively to clean greasy and dirty hands also it kills about 99.9% of microbes, as well as, it contains glycerin & Vitamin E for moistening the skin. Liquid Florosept-handwash makes skin soft and very clean with pleasant odor.

### :وصف المنتج

منظف الايدي فلوروسبت هو منظف جاهز للاستخدام له خصائص تنظيف عالية الكفاءة ومضاد للجراثيم. يتميز بكفائته العالية مع تنظيف الايدي المتسخة وال لترطيب الجلد وجعله ناعماً. منظف الاي دي من فلوروسبت يترك الايدي نظيفة جيدا وناعمة مع رائحة مميزة

### **Active Ingredients**

Aqua, auaternary ammonium compound, Sodium lauryl sulphate, Sodium chloride, C ocamidepropyl Betaine, Cocomide DEA, Glycerin, Vitamin E, Parfum, tetra sodium EDTA, citric acid, color, perfume, preservative & auxiliary materials.



1000 ml

ASC-FLSP-HS-B - 1000ml

# FLÓRÓSEPT

**ANTIBACTERIAL** 

### Hand Disinfectant

Chlorhexidine gluconate 2 % W/V.
Isopropyl Alcohol 4 %, W/V

### **Product Description**

Liquid Florosept-handwash is a ready to use, detergent-based hand wash that provides superior cleaning and anti-bacterial effects. Florosept-handwash is strong and very effectively to clean greasy and dirty hands also it kills about 99.9% of microbes, as well as, it contains glycerin & Vitamin E for moistening the skin. Liquid Florosept-handwash makes skin soft and very clean with pleasant odor.

### :وصف المنتج

منظف الايدي فلوروسبت هو منظف جاهز للاستخدام له خصائص تنظيف عالية الكفاءة ومضاد للجراثيم، يتميز بكفائته العالية مع تنظيف الايدي المتسخة وال لترطيب الجلد وجعله ناعماً، منظف الاي دي من فلوروسبت يترك الايدي نظيفة جيدا وناعمة مع اللحة ممنة

### Active Ingredients -

Aqua, quaternary ammonium compound, Sodium lauryl sulphate, Sodium chloride, C ocamidepropyl Betaine, Cocomide DEA, Glycerin, Vitamin E, Parfum, tetra sodium EDTA, citric acid, color, perfume, preservative & auxiliary materials.

1000 ml

ASC-FLSP-HS - 1000ml





# ASC FABRIC SOFTENER

منعم ومعطر للملابس

### **Product Description**

ASC- fabric softener Offer (give) your clothes an amazing freshness and softness with ASC- fabric softener, makes ironing easier, with long-lasting fragrance.

:وصف المنتج

يعطي أي أس ُسي- منعم الاقمشة لملابسك انتعاشا و نعومة فائقة , يجعل الكي اكثر سهولة , مع رائحةدوم ثيرا

### Ingredients

Cationic surfactants dipalmitoylethyl hydroxyethylmonium methosulfate, nonionic surfactants, caprylyl glucosid, aqua, glycerin, camellia sinensis leaf extract, ethylene brassylate, ethyl linalool, methyl dihydrojasmonate, benzyl acetate, dimethyl heptenal, fragrance

30 Litter	ASC-FS-30L		
4 Litter	ASC-FS-4L		
4 Lillei	A3C-1 3-4L		



# ASC GLASS CLEANER

منظف وملمع زجاج بتركيبة فعالة ومتطورة

### **Product Description**

Easily removes sediment, dirt and fingerprints Gives crystal transparency and streak-free shine ldeal for cleaning glass and other polished surfaces Expiration, formulated to clean quickly without ammonia. Made without parabens, phosphates and phthalates Readily biodegradable, pH balanced.

وصف المنتج

يزيل الرواسب والأوساخ وبصمات الأصابع بسهولة يمنح ش افية كريستالية ◘₪ خا من الخطوط ل ف نً مثالي لتنظيف الزجاج والأسطح الأخرى , صمم للتنظيف السريع , خا من الأمونيا ، البرابين ، ل ال و سات وال يثال ت. قا ل للتح ل البيولو ي ب هولة ، در ة الحمو ة متعاد ة

### **Active Ingredients -**

alcohol, surfactant, blue color, fragrance, water, Antistatic additives

1000 ml

ASC-GC - 4L





# ASC-POVIDONE IODINE 10 % W/V

أي أس سى- بوفيدون أيودين ١٠ ٪ و/ح

### Product Description -

ASC-Povidone Iodine %10 Antiseptic Solution is a broad-spectrum bactericidal thatis effective against gram (+) and gram (-) bacteria, viruses, fungi & spores. The effect beginning with the film formation on the wound, continues for a long time by the releasing of controlled iodine from the film layer. ASC-Povidone Iodine %10 Antiseptic Solution is used as a skin and mucosa antiseptic before surgery, injection, biopsy, blood collection and ponction operations and for wound treatment, pyoderma, dermatosis and as hygienic and surgical hand antiseptic.

### وصف المنتج

### · Product Composition -

Povidone iodine %10 w/v.

### Product Indication

Shows effective microbicidal activity against both gram-positive and gram-negative bacteria, viruses, fungi, protozoa, spores, rickettsia and yeasts.

250 ml

ASC-PVP-801-250 ml



### **Product Description** -

 ${\sf ASC-OPA/28} \ is \ an \ or tho-Phthalal dehyde \ (OPA) \ Solution \ High-Level \ Disinfectant \ for semi-critical \ medical \ devices.$ 

### **Active Ingredients**

**Active Ingredient** 

ortho-Phthalaldehyde (OPA)... 0.575%

Inactive Ingredients...... 99.425%

**Buffers** 

Surfactants

Antifoaming agent

Solvent

Deionized water

Total......100.0%

### **Effective Spectrum**

Destroy pathogens in just 10 minutes (at room temperature) with this powerful, ready-to-use disinfectant. Some reprocessors of semi-critical medical devices use the active ingredient, 0.575% ortho-phthalaldehyde (OPA), to safely disable harmful pathogens, like:2

- MRSA
- VRE
- TB
- HIV
- Hepatitis Viruses



5Litter | 50 strips

ASC-OPA-005S



# ASC-CLOREX Shower Bath Kit

Pre-operative antiseptic shower bath-kit with Chlorhexidine Gluconate (CHG) Antiseptic solution

### Product Description

ASC-Shower bath-KIT is a new product line developed for hospitals to reduce the risk of surgical site infections. Bathing with CHG has been shown to reduce bacterial level, which can help to prevent of surgical site from occurring infections and thus help to reduce surgical complications. CHG also has excellent antiseptic properties and has been proven to be superior to other preparation products for bacterial decolonization. CHG rapidly kill transient and resident microorganisms on the skin and is effective against Gram-Positive and Gram-Negative bacteria.

### **Ingredients**

Composition 100g: 4.0g, coformulants and water up to 100g. chlorhexidine gluconate (CHG)

KIT

ASC-SHB-201-100 ml



# FLÓRENE ANTIBACTERIAL DEODORANT BAR SOAP

DEEP CLEANING SOAP THAT WASHES AWAY HARMFUL GERMS & BACTERIA

### Product Description -

FLÓRENE ANTIBACTERIAL DEODORANT BAR SOAP is designed to provide excellent germ protection. Contains antibacterial agents that eliminate bacteria and moisturizes that improve the skin's texture. Deodorant leaves the skin smelling fresh and rids the body of bad odors.

وصف المنتج: تم تصميم فلورين صابون مضاد للبكتريا لتوفير حماية ممتازة من الجراثيم. يحتوي على عوامل مضادة للجراثيم تقضي على البكتيريا وترطب وتحسن نسيج الجلد. مزيل للعرق يترك البشرة برائحة منعشة ويخلص الجسم من الروائح الكريهة.

### Active Ingredients -

FLÓRENE Active Ingredient - Purpose: BENZALKONIUM CHLORIDE 0.10% - Antibacterial. Inactive Ingredients: Soap (Sodium Cocoate\*, Sodium Palm Kernelate\*, Sodium Palmate\*), Water, Talc, Fragrance, Glycerin, Sorbitol & Auxiliary materials.



30g

ASC-ADBS-30g



### **ACHIEVEMENT SCOPE COMPANY FOR MEDICAL INDUSTRIES** شركة نطاق الإنجاز للصناعات الطبية



المملكة العربية السعودية KINGDOM OF SAUDI ARABIA







